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Dear customer,

Our first aim is to provide a medical device in constant high quality. During our permanent survey and control in the field, we have observed an error that could prevent a reliable re-start of the **corpuls³** device used in your company.

We kindly request you to pay attention to the warning for an unintended turn-off of the equipment as well.

Please read this customer information carefully and return the filled in and duly signed confirmation letter (Annex B1) to GS until **25.10.2007** at the latest.

The serial numbers of all affected devices are mentioned in Annex C. The serial numbers of the devices used in your company – according to our company documentation – are mentioned in Annex D.

The error description as mentioned under both items is only limited to those **corpuls³** devices that have been supplied until delivery date 19.10.2007. All other devices of the manufacturer GS Elektromedizinische Geräte G. Stemple GmbH are not affected by this safety information.

1. Safety information - Error prevents restart of device

1.1 Error description

In case the **corpuls³** operator presses the green ON/OFF button during the power-off sequence (or directly after) of the compact device once again, the Monitoring Unit starts without the Patient Box and/or the Defibrillator/Pacer Unit. Hence the compact device is not completely turned on and not ready for operation.

On the screen of the Monitoring Unit this situation is displayed as an irregular connection by 2 flashing bars between the symbols of the modules.


Fig. 1: Indicator during failure appearance

1.2 Troubleshooting

In case of the failure indication as mentioned under item 1.1 (see Fig. 1), the failure can be safely rectified by a temporary separation and re-connection of the modules, as the turned off modules can now be powered up again under control. Repeated turning off of the Monitoring Unit is not necessary.

1.4 Immediate measure

We kindly request you to inform all **corpuls[®]** operators within your organisation at short notice NOT to press the green ON/OFF button of the Monitoring Unit (dark screen) at least 3 seconds after the power-off sequence. Kindly refer the operator to the corrective action after appearance to this failure.

1.4 Manufacturer action

The mentioned error is fixed in a new firmware, available as from 11/2007. Implementation of the firmware update is carried out by the manufacturer (GS Elektromedizinische Geräte G. Stemple GmbH) or its authorised service agents and distributors.

If no service contract is signed with one of our authorised **corpuls[®]** service agents or distributors and the maintenance of the equipment is carried out by third parties, don't fail to contact the authorised **corpuls[®]** service agent or distributor immediately in order to make an appointment for this important firmware update.

1.5 Deadline

The implementation of this corrective action has to take place until **01.05.2008** at the latest. For all **corpuls³** devices without this updated firmware the manufacturers guarantee expires on this date and an operation is no longer allowed.

2. Safety information - Unintentional shutdown of the device

2.1 Error description

In the course of our research we have noticed that during operation of a **corpuls³** in few cases an accidental and unnoticed shutdown caused by an unintentional operating error is carried out. In doing so the user has unintentionally pressed the green ON/OFF button and did not notice and interrupt the consequential countdown.

1.4 Immediate measure

We kindly request you to inform all users within your organisation on short term, that in case of unintentional operation of the green ON/OFF button, a dialog box with countdown appears on the Monitoring Unit. If the operator does not interrupt the countdown within 10 seconds, the unit is turned off and can only be turned on after the complete power-down (see 1.1).

2.3 Manufacturer action

As from 11/2007 the software version 1.3.0 is going to be available which also includes a new shutdown procedure. The update is carried out together with the firmware update by the manufacturer or its authorised service agents or distributors. As from this time the shutdown of the device must be confirmed by operation of a second button.

If no service contract is signed with one of our authorised **corpuls[®]** service agents or distributors and the maintenance of the equipment is carried out by third parties, don't fail to contact the authorised **corpuls[®]** service agent or distributor immediately in order to make an appointment for this important software update.

2.4 Deadline

The implementation of this corrective action has to take place until **01.05.2008** at the latest.

We count on your understanding for the realisation of this quality assuring procedure.

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Contact person for inquiry calls: For service agent in charge or the local distributor see <http://www.corpuls.com>.

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Carsten Fuchs
Vice President, Customer Support
Serviceleiter

Markus Raab
Quality Management
Qualitätsmanagement

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


Annex B 1

- reply form GB -

Please cross all fields which apply your company.

- Hereby we confirm to be read and understood the safety information of GS Elektromedizinische Geräte G. Stemple GmbH of 2007-10-19.
- We have following **corpuls³** devices:

No	Monitoring unit	Patient box	Defibrillation / Pacer unit
1			
2			
3			
4			
5			

- We have briefed the users in our company / organisation about the shutdown sequence and the required delay before restart.
- We run out of **corpuls³**, whose serial number is listed in Annex C of this document. (Please describe what was done with these devices, e.g: sort out, scrapped, sold to If possible, please attach a copy of the cross-reference.)

Filled in by customer (please in block letters):

Company / Organisation: _____

Address: _____

City: _____

Country: _____

Name: _____

First name: _____

Title: _____

Fax: _____

Phone: _____

Company stamp: _____

e-mail-adress: _____

Date / Signature: _____

Please fill in this reply form and send or fax it to us by **25.10.2007**.

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Hans Peter Esser GmbH
Meßmer Medizintechnik GmbH