

**Safety Notice**  
**Technical Bulletin no. 005**



GS Elektromedizinische Geräte  
G. Stemple GmbH  
Hauswiesenstraße 26  
D-86916 Kaufering  
Tel. +49 8191 65722-0  
Fax +49 8191 65722-22  
info@corpuls.com  
www.corpuls.com

Nr. 005	Target audience all users	Date 2008-12-12	Nr of pages 8
Concerned products <b>corpuls<sup>3</sup></b>	Serial numbers / Lot identification all devices	Software / Firmware each software version up to and including 1.4.2	

Dear user,

With this letter we are informing you of a **corpuls<sup>3</sup>** Defibrillator/Monitoring system series software issue.

This letter is intended to inform you about

- what the problem is and under what conditions it may appear,
- the measures, you as operator must take to minimize the effects of this problem,
- the measures scheduled by the manufacturer to sustainably fix this problem.

Affected are all **corpuls<sup>3</sup>** devices independent of both configuration and software version up to and including software version 1.4.2. We kindly request you to read this Safety Notice carefully and return the completed and signed Confirmation Letter annex B not later than January 31, 2009.

## 1. Description of error



Under appearance of the below mentioned conditions the software may store data in a non-provided memory area. So far, thereupon the following malfunctions have been noticed:

- The message "Pairing failed" may appear after start pairing of the different device components Defibrillator Unit, Patient Box and Monitoring Unit.
- The printer may fail when the speed is changed.
- Temperature measurement may suddenly fail.

## 2. Precondition for the appearance of the error

The error may appear under simultaneous appearance of the below listed conditions:

1. The CompactFlash<sup>®</sup> memory card is full, almost full, damaged or being removed during operation.
2. Several diagnosis ECGs are being recorded.
3. Subsequently, configuration data (e. g. pairing data, alarm limits) are being modified.

Document name and location:	U:\Allgemein\Technische Bulletins\TB_005_EN.doc	Page	1 of 8	 
Creation date:	2008-12-11	Release date:	2008-12-15	
Originator name:	Carsten Fuchs	Release name:	Klaus Stemple	

### 3. Potential risk

The Monitoring Unit indicates the message „Pairing failed“. The connection between the device components Defibrillator Unit, Patient Box and Monitoring Unit is interrupted. The device is not ready for use.

### 4. Safety information

Until installation of a new software version please consider the items below:

- Archive data periodically from the CompactFlash® memory card, than delete them completely; latest if the message “CF card capacity low” appears (see User Manual, page 203).
- Do not remove the CompactFlash® memory card from the device during operation.
- Do not generate further diagnosis ECGs after any error message concerning the CompactFlash® memory card has appeared during operation. Contact contact the technical customer service of an authorized sales and service partner.
- Do not change any configuration adjustments (User Manual chapter 7) during regular operation.
- Do not start a pairing process during regular operation.
- Perform the functional test according to the User Manual (chapter 9.2.1) periodically.

### 5. Trouble shooting at devices which caught attention

If the aforementioned error has already been noticed, we recommend the reinstallation of the software currently operated in your device until the advertised software version 1.4.3 is available. Contact the technical customer service of your authorised sales and service partner.

### 6. Immediate measures



Please instruct all users in your organisation without delay about the potential risk and measures. If doubts exist, the measures described in the User Manual, chapter 10.2 “Troubleshooting and remedying malfunctions” must be executed immediately.

### 7. Manufacturers measures

Starting with calendar week 02/2009 a new software version will be provided. The version number is 1.4.3.

The software installation is carried out by the manufacturer or by an authorised sales and service partner.

If you don't have any service contract with an authorised **corpuls**® sales and service partner, and your device is serviced by third party companies, please do not hesitate to immediately contact your national sales and service partner to make an appointment for software installation.

Document name and location:	U:\Allgemein\Technische Bulletins\TB_005_EN.doc	Page	2 of 8		
Creation date:	2008-12-11	Release date:	2008-12-15		
Originator name:	Carsten Fuchs	Release name:	Klaus Stemple		



**Safety Notice**  
**Technical Bulletin no. 005**



**8. Deadline**

The measures must be implemented not later than May 01, 2009. All **corpuls<sup>3</sup>** devices with software version 1.4.2 or lower (indicated in the menu System – Info) are not allowed for operation as from this date.

System - Info		
System	Software Rev.	Seriennummer
Monitoreinheit	REL-1.3.0_C3_BP	A10000002BF35D05
Patientenbox	REL-1.3.0_C3_BP	7D0000002BFFF505
Defibrillator	REL-1.3.0_C3_BP	190000002C15F205
Optionen	Software Rev.	Seriennummer
Biphas. Modul	M:v2.00K/S:v2.00G	--
EKG	1M	--
SpO2	1/v4.5.0.3/v1.0.0.4	--
NIBD	LM3.390/SM V220/0533	--
CO2	--	--
IBD	v10/v5/v5	--
Temp	v8/v4/v4	--
GSM	MC55 04.00	--
EKG-Interpretation	18.24-03	--
<b>Info</b>		Abbr. OK

*Indication of software version*

May we thank you for your understanding and apologise for any inconvenience you may have in connection with this corrective action. Questions concerning this matter will be answered by your national sales and service partner (see annex D or [www.corpuls.com](http://www.corpuls.com)).

Contact person of the manufacturer for inquiry calls:


Mr. Carsten Fuchs,  
 Vice President, Customer Support



Phone: +49 81 91 6 57 22 30  
 Fax: +49 81 91 6 57 22 22  
 E-Mail: [fuchs@corpuls.com](mailto:fuchs@corpuls.com)

With kind regards  
 GS Elektromedizinische Geräte G. Stemple GmbH

  
 Günter Stemple  
 Managing Director  
 Geschäftsführer

  
 Klaus Stemple  
 General Manager R&D/Production  
 Geschäftsführer F&E/Fertigung


  
 Carsten Fuchs  
 Vice President, Customer Support  
 Serviceleiter

Document name and location:	U:\Allgemein\Technische Bulletins\TB_005_EN.doc	Page	3 of 8	 
Creation date:	2008-12-11	Release date:	2008-12-15	
Originator name:	Carsten Fuchs	Release name:	Klaus Stemple	

**Annex A**

Figure of device combination **corpuls<sup>3</sup>**

- 1 – Monitoring Unit
- 2 – Patient Box
- 3 – Defibrillator



*Type plates with positions of serial numbers*

**Safety Notice**  
**Technical Bulletin no. 005**



**Annex B**

Reply form

Kindly mark all fields where applicable for your organization.

- We have read and understood the Safety Notice of GS Elektromedizinische Geraete G. Stemple GmbH from December 12, 2008.
- We operate the below listed devices of the **corpuls<sup>3</sup>** series (pls. enter all serial numbers):

No.	Monitoring Unit	Patient Box	Defibrillator Unit
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			

- All our employees have been conveniently informed about the content of this Safety Notice.
- We have discontinued the operation of the following **corpuls<sup>3</sup>** whose serial number is listed in annex C of this document. Kindly specify what was done with this device; e.g. decommissioning, selling, other use. If available, please attach a cross reference.

---



---



---



---

To be completed by the customer (block letters please):

Organization: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

Country: \_\_\_\_\_

Name: \_\_\_\_\_

First name: \_\_\_\_\_

Title: \_\_\_\_\_

Fax: \_\_\_\_\_

Phone: \_\_\_\_\_

Company stamp: \_\_\_\_\_

e-mail-address: \_\_\_\_\_

Date / Signature: \_\_\_\_\_

Please complete this reply form and return by fax not later than **31-01-2009**.

GS Elektromedizinische Geraete G. Stemple GmbH, Hauswiesenstrasse 26, D-86916 Kaufering

**Fax: + 49 8191 65722 - 22**

Document name and location:	U:\Allgemein\Technische Bulletins\TB_005_EN.doc	Page	5 of 8	
Creation date:	2008-12-11	Release date:	2008-12-15	
Originator name:	Carsten Fuchs	Release name:	Klaus Stemple	



## Annex C

### National sales and service partners

#### Czech Republic

CHEIRON a.s.  
Blatenska 27 a  
CZ 32600 Plzen  
phone: +420-377 590 412  
fax: +420-377 590 435  
e-mail: [vmattasova@ceiron.cz](mailto:vmattasova@ceiron.cz)

#### Estonia

saarik & ko  
Akadeemia tee 33  
EE 12618 Tallinn  
phone: +372 6525646  
fax: +372 6525388  
e-mail: [saarik@saarik.ee](mailto:saarik@saarik.ee)

#### Hungary

ANAMED  
Analytical Medical Instruments Kft.  
Köszeg u. 29.  
HU 1144 Budapest  
phone: +36 1 2209236  
fax: +36 1 2215531  
e-mail: [fischer@anamed.hu](mailto:fischer@anamed.hu)

#### Iran

Eshtood Kar Co. Ltd.  
Unit 7, No. 4 Golshahr Blvd. Africa Ave.  
IR11369 Tehran  
phone: +98-21 22040199  
fax: +98-21 22 053481  
e-mail: [eshtood@yahoo.com](mailto:eshtood@yahoo.com)

#### Israel

ARDON Mediacal Equipment Ltd.  
24 Ha'Charoshet Street (PO Box 378)  
Or Yehuda 60375  
phone: +972 3 5333236  
fax: +972 3 5334801  
e-mail: [arie@ardon.co.il](mailto:arie@ardon.co.il)

#### Italy

Mortara Rangoni EUROPE s.r.l  
Via Cimarosa, 103/105  
40033 Casalecchio di Reno (BO)  
phone: +39 051 2987811  
fax: +39 051 6133582  
e-mail: [rabito@mortara.it](mailto:rabito@mortara.it)

#### Kingdom of Saudi Arabia

Scientific & Medical Equipment House  
P.O. Box 15 84  
11441 Riyadh  
phone: +966 1 4647711  
fax: +966 1 4631507  
e-mail: [jerome@smeh.com.sa](mailto:jerome@smeh.com.sa)

#### Netherlands, Belgium and Luxembourg

corpuls Nederland B.V.  
Chr. Huygensweg 25A  
NL 3225 LD Hellevoetsluis  
phone: +31 181 390963  
fax: +31 181 390970  
e-mail: [kkamphuis@corpuls.nl](mailto:kkamphuis@corpuls.nl)

#### Poland

TEHAND Ltd  
ul. Przybyszewskiego 17B  
PL 30-128 Kraków  
phone: +48 126610180  
fax: +48 126620354  
e-mail: [zwolakm@tehand.pl](mailto:zwolakm@tehand.pl)

#### Romania



Deltamed Ltd.  
Calea Manastur 70 / 30  
RO-400658 Cluj-Napoca, Jud. Cluj  
phone: +40 264 427666  
fax: +40 264 427679  
e-mail: [dan.gorgan@dgis.ro](mailto:dan.gorgan@dgis.ro)

#### Spain

Diagniscan - Critical Care Division  
CH Werfen - Hospital Group  
Aragón, 90 ; 08015 Barcelona  
phone: +34 93 401 0376  
fax: +34 93 401 0154  
e-mail: [jrnavarro@diagniscan.es](mailto:jrnavarro@diagniscan.es)

#### Sultanate of Oman

Seven Seas Co. LLC.  
P.O. Box 1222  
C.P.O Seeb-121  
phone: +968-24533471  
fax: +968-24533472  
e-mail: [medical@sevenseasoman.com](mailto:medical@sevenseasoman.com)

Document name and location:	U:\Allgemein\Technische Bulletins\TB_005_EN.doc	Page	6 of 8	 
Creation date:	2008-12-11	Release date:	2008-12-15	
Originator name:	Carsten Fuchs	Release name:	Klaus Stemple	

**Safety Notice**  
**Technical Bulletin no. 005**



**Switzerland**

Schiller-Reomed AG  
Riedstraße 14  
CH 8953 DIETIKON  
phone: +41 44 7410209  
fax: +41 44 7403710  
e-mail: [ernst@schiller-reomed.ch](mailto:ernst@schiller-reomed.ch)

**Turkey**



EMS Mobile Systems and  
Hospital Equipments Industry & Trade Inc.  
1.Organize Sanayi Bölgesi  
Uygurlar Cad. No: 5/A  
TR 06930 SINCAN-ANKARA  
phone: +90 212 511 04 34  
fax: +90 212 527 61 11  
e-mail: [mozgun@ems.tc](mailto:mozgun@ems.tc)

**Other territories**

Weinmann Geräte für Medizin  
GmbH & Co. KG  
Kronsaalsweg 40  
22525 Hamburg  
phone: +49 40/54702-337  
fax: +49 40/54702-467  
e-mail: [m.szepannek@weinmann.de](mailto:m.szepannek@weinmann.de)

**Manufacturer:**

GS Elektromedizinische Geraete G. Stemple GmbH  
Hauswiesenstrasse 26  
D-86916 Kaufering  
phone: +49 8191 65722-0  
fax: +49 8191 65722-22  
e-mail: [info@corpuls.com](mailto:info@corpuls.com)  
[www.corpuls.com](http://www.corpuls.com)

Document name and location:	U:\Allgemein\Technische Bulletins\TB_005_EN.doc	Page	7 of 8	 
Creation date:	2008-12-11	Release date:	2008-12-15	
Originator name:	Carsten Fuchs	Release name:	Klaus Stemple	

**Safety Notice**  
**Technical Bulletin no. 005**



**Annex D**

Serial numbers

Ladies and Gentlemen,

According our records you're operating the **corpuls<sup>3</sup>** device with the serial numbers as listed below (for position of serial number see annex A):

Monitoring Unit	Patient Box	Defibrillator		Monitoring Unit	Patient Box	Defibrillator